Leading Edge Custom Queries and their Applications Across Ontario

Hope Latam, QIDS Specialist

AFHTO Conference

Toronto, ON

Oct 28th – 29th 2015



Presenter Disclosure

- Presenters: Hope Latam
- Relationships with commercial interests:
 - Not Applicable

Disclosure of Commercial Support

• No Commercial Support

- Potential for conflict(s) of interest:
 - None

Mitigating Potential Bias

Not Applicable

Custom Forms


Entering data in PS:

- Free text
- Stamps
- Encounter Assistants
- Custom Forms






Structured and organized data input




Extracting Data from Custom Forms

Basic searches don't work...



Tried exporting search results using flowsheets...



Custom Queries

Edit Custom Queries

FROM STRUCTURED STAMP ITEM VALUE iv

sv.structured stamp version id


as atmostured atoms meater id

Custom Form Named Items


SELECT

Worked with TELUS and Cognisant MD to develop a set of custom queries

JOIN STRUCTURED STAMP ss ON sv.structured stamp master id =



Deployment



Edit Custom Queries

Custom Form Named Items

SELECT

iv.patient_piece_id, pp.patient_id, pp.patient_piece_date, si.item_name
iv.stamp_item_key, DBMS_LOB.substr(iv.stamp_item_value,256) AS ItemValue
FROM STRUCTURED_STAMP_ITEM_VALUE iv

JOIN PATIENT_PIECE pp ON pp.patient_piece_id = iv.patient_piece_id


JOIN STRUCTURED_STAMP_VERSION sv ON pp.structured_stamp_note_stamp_id :
sv.structured stamp_version_id

JOIN STRUCTURED_STAMP_ITEM si ON si.stamp_item_key = iv.stamp_item_key si.structured_stamp_version_id = sv.structured_stamp_version_id JOIN STRUCTURED_STAMP ss ON sv.structured_stamp_master_id =

as atmostured atoms meater id




Custom Query Converter



	Α	В	С	D	Е	F	G
1	PATIENT_PIECE_ID 💌	PATIENT_ID 🔻	PATIENT_PIECE_DATE	Smoker	Ex Smoker 🔻	Referred T	Quit Date
2	Form 1	123	1-Apr-15	TRUE	FALSE	TRUE	1-May-15
3	Form 2	345	20-Apr-15	FALSE	TRUE	FALSE	1-Jan-14
4							
5							
6							
7							


Example




PS Suite® EMR (1387-Drayton Community Health Centre)

File Settings Reports Messages Window Help

Example





	Α	В	С	D	E	F	G	Н	1	J
1	PATIENT_PIECE_ID 💌	PATIENT_ID 🔻	PATIENT_PIECE_DATE	CurrentSmoker 💌	ExSmoker 💌	NeverSmoked ▼	WantQuitYes 🔻	WantQuitNo 🔻	WantQuitNA 🔻	ReferToCounsellorYes 🔻
2	1534419	20802	2-Mar-15	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE
3	1534767	6815	2-Mar-15	FALSE	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE
4	1534946	26075	2-Mar-15	FALSE	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE
5	1535228	11417	2-Mar-15	FALSE	FALSE	TRUE	FALSE	FALSE	FALSE	FALSE
6	1535859	1803	3-Mar-15	FALSE	FALSE	TRUE	FALSE	FALSE	FALSE	FALSE
7	1535882	4857	3-Mar-15	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE
8	1536030	17002	3-Mar-15	FALSE	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE
9	1536381	15173	3-Mar-15	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE
10	1536382	28844	3-Mar-15	FALSE	FALSE	TRUE	FALSE	FALSE	FALSE	FALSE
11	1536601	8385	3-Mar-15	FALSE	FALSE	TRUE	FALSE	FALSE	FALSE	FALSE
12	1536645	27951	3-Mar-15	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE
13	1536777	21605	3-Mar-15	FALSE	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE
14	1536787	20477	3-Mar-15	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE	FALSE

Example

	А	В	С	D	Е	F	G	Н	1	J
1	PATIENT_PIECE_ID 🔻	PATIENT_ID 🔻	PATIENT_PIECE_DATE	CurrentSmoker	▼ ExSmoker ▼	NeverSmoked ▼	WantQuitYes 🔻	WantQuitNo 🔻	WantQuitNA 🔻	ReferToCounsellorYes 🗷
7	1535882	4857	3-Mar-15	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE	TRUE
9	1536381	15173	3-Mar-15	TRUE	FALSE	FALSE	FALSE	FALSE	FALSE	TRUE
17	1537271	15173	4-Mar-15	TRUE	FALSE	FALSE	TRUE	FALSE	FALSE	TRUE
28	1544923	26374	11-Mar-15	TRUE	FALSE	FALSE	TRUE	FALSE	FALSE	TRUE
59	1552341	27337	19-Mar-15	TRUE	FALSE	FALSE	TRUE	FALSE	FALSE	TRUE
82	1555856	20207	24-Mar-15	TRUE	FALSE	FALSE	TRUE	FALSE	FALSE	TRUE
3074										
3075										
3076										,
3077										
3078										,
3079										,
3080										,
3081										
3082										
3083										
3084										
3085										
3086										
3087										
3088										
3089										
3090										
3091										
3092										

Success Stories

Prior to the Custom Query, I would spend hours creating searches in the EMR. This was very time consuming and I was not always confident in the accuracy of these searches. With the Custom Query, I can enter the name of the custom form and exact date range I want to search. Within seconds I have all of the data fields extracted and ready to be analyzed. The Custom Query Converter is like the cherry on top! It takes out all the manual work of sorting the columns, which allows for the data to be analyzed even faster. All of our custom form searching needs can be accomplished using this handy tool. I am very grateful to be able to have the Custom Query in our EMR and be able to provide my team with true, relevant data.

-Desiree Leslie, New Vision FHT

Success Stories

Thank you EWFHT for developing the custom queries and the excel converter! We have just recently completed our quarter 4 IHP stats and used these new tools quite extensively. Previously, it would take the better part of a week to tabulate quarterly IHP stats but with the help of these queries and some custom form development we have managed to reduce the time to just a few hours. AMAZING!

-Joe Szeplaki, Garden City FHT

